

**Nome Prodotto**















Carrozzina ad Autospinta Viky Light.

**Destinazione d'uso**

Le carrozzine sono un mezzo di trasporto indispensabile per le persone disabili con assente o ridotta capacità di deambulazione.

Lo scopo è quello di facilitare gli spostamenti autonomi e di salvaguardare il benessere fisico della persona, evitando la comparsa di complicanze correlate alla posizione seduta obbligatoria (piaghe da decubito, dolori posturali cronici...).



Codice	Taglia	Lunghezza	Larghezza	Altezza	Larghezza Piegata	Larghezza Seduta	Profondità Seduta	Altezza Schienale	Altezza Seduta da Terra	Portata Massima (Kg)	Dimensioni Ruote Anteriori	Dimensioni Ruote Posteriori	Peso Netto con Appoggiatesta	Dimensioni Scatola
Misure in Cm														
01.04.630DQR41FD	41	80	61	90.5	29.5	41	40	44	50	120	20 (8")	60 (24")	13	82x31x94
01.04.630DQR43FD	43	80	63	90.5	30	43	40	44	50	120	20 (8")	60 (24")	13.1	82x31x94
01.04.630DQR46FD	46	80	65	90.5	30.5	46	40	44	50	120	20 (8")	60 (24")	13.2	82x31x94
01.04.630DQR48FD	48	80	67	90.5	31	48	40	44	50	120	20 (8")	60 (24")	13.2	82x31x94

**Caratteristiche**

Telaio	Alluminio Verniciato a Polvere
Crociera	Doppia
Colore	Grigio Canna di Fucile (Opaco)
Ruote Anteriori	20 x 2,5 cm (8" x 1") PU Grigie
Ruote Posteriori	Estraibili da 60 cm (24") PU Grigie
Anello Corrimano	Alluminio
Accessori	Ruote Anti-ribaltamento e Ruotine per Ascensore
Braccioli	Ribaltabili, Modello a Scrivania con Poggiabraccio in PU
Pedane	Estraibili e Regolabili con Supporto per Gambe
Poggiapiedi	PVC Nero
Sedile e Schienale	Nylon Morbido Imbottito, lavabile, disinfettabile e ignifugo

**Manutenzione e Pulizia**

Disinfezione possibile con disinfettanti diluiti al 5% anche a base alcolica

**Conservazione e Smaltimento**

Si consiglia l'uso ad un solo utente

Conservare al riparo da fiamme libere e fonti di calore

Smaltire il dispositivo rispettando le normative vigenti in materia di tutela dell'ambiente e raccolta differenziata

**Confezione e Nr. Articoli**

Scatola	1 Pezzo
---------	---------

**Garanzia**

Europea	24 mesi per difetti di fabbricazione ad esclusione delle parti consumabili
---------	--

**Riferimenti del Prodotto**

Fabbricante	Demarta Virginio S.r.l.
Classe di Dispositivo MDR 2017/745	Classe I
Codice ISO DM332/99	12.21.06.033 / 12.21.06.039 / 12.24.06.175
Codice ISO DPCM 12.01.2017	12.22.03.006
Codice CND	Y122106

Codice RDM

1567089/R















### Product Name

Manual Wheelchairs Viky Light.

### Use

The wheelchairs are an essential way of transport for disabled people with none or limited walking ability. The purpose is to facilitate autonomous movements and to safeguard the physical well-being of the person, avoiding the appearance of complications related to the sitting position (bed sores, chronic postural pains ...).



Code	Size	Length	Width	Height	Width Folded	Seat Width	Seat Depth	Back Height	Floor to Seat Height	Weight Capacity	Front Wheel	Rear Wheel	Weight with Footrest	Carton Size
Measures in Cm														
01.04.630DQR41FD	41	80	61	90.5	29.5	41	40	44	50	120	20 (8")	60 (24")	13	82x31x94
01.04.630DQR43FD	43	80	63	90.5	30	43	40	44	50	120	20 (8")	60 (24")	13.1	82x31x94
01.04.630DQR46FD	46	80	65	90.5	30.5	46	40	44	50	120	20 (8")	60 (24")	13.2	82x31x94
01.04.630DQR48FD	48	80	67	90.5	31	48	40	44	50	120	20 (8")	60 (24")	13.2	82x31x94

### Features

Frame	Powder Coating Aluminum Frame
Cross Bar	Double Cross Brace
Color	Smokey Grey (Dull Polish)
Front Wheels	20 x 2,5 cm (8" x 1") PU Grey Tyre
Rear Wheels	60 cm (24"), Spoke Wheels, PU Grey Tyre, Quick Release
Handrim	Aluminum
Accessory	Small Wheels for Lift and Anti Tipper Wheels
Armrest	Flip Back Slope Model Armrest with PU Armpad
Footrest	Swing In and Away, Length Adjustable with Leg Support
Footplate	Plastic Black
Seat and Backrest	Soft and Padded Nylon, washable, disinfectable and fire-retardant

### Maintenance and Cleaning

Disinfection is possible with 5% diluted disinfectants even with alcohol base

### Storage and Disposal

Single-user Product

Keep away from heat, sparks, open flame or other sources of ignition

Dispose of the device in compliance with current regulations on environmental protection and separate waste collection

### Package and Nr. Items

Box 1 Piece

### Warranty

European	24 months for manufacturing defects with the exception of consumable parts
<b>Product References</b>	
Manufacturer	Demarta Virginio S.r.l.
Directive MDR 2017/745	Class I
ISO Code DM332/99	12.21.06.033 / 12.21.06.039 / 12.24.06.175
ISO Code DPCM 12.01.2017	12.22.03.006
CND Code	Y122106
RDM Code	1567089/R