

**Nome Prodotto**

Protesi Maxima Simmetrica Cava.



Codice	Taglia	Taglia Francese	Taglia Europea
06.01.4001C.0	0	65	50
06.01.4001C.1	1	70	55
06.01.4001C.2	2	75	60
06.01.4001C.3	3	80	65
06.01.4001C.4	4	85	70
06.01.4001C.5	5	90	75
06.01.4001C.6	6	95	80
06.01.4001C.7	7	100	85
06.01.4001C.8	8	105	90
06.01.4001C.9	9	110	95
06.01.4001C.10	10	115	100
06.01.4001C.11	11	120	105
06.01.4001C.12	12	125	110
06.01.4001C.13	13	130	115
06.01.4001C.14	14	135	120

Destinazione d'uso

Le protesi mammarie esterne vengono utilizzate in soggetti che hanno subito una mastectomia. Grazie a questi dispositivi la morfologia del corpo resta invariata.

Descrizione

La PROTESI MAXIMA SIMMETRICA CAVA è composta da un gel di silicone avvolto da un sottile film di poliuretano; concepita per ripristinare una visione anatomico-estetica del corpo. La protesi è anallergica ed esteticamente naturale, la forma a goccia la rende adatta ad ogni tipo di seno.

Materiali	
Protesi	Silicone
Manutenzione e Pulizia	
La pulizia è possibile con un panno umido e un detergente neutro	
Si consiglia di utilizzare deodoranti, profumi o altri cosmetici prima di indossare le protesi	
Conservazione e Smaltimento	
Si consiglia l'uso ad un solo utente	
Conservare al riparo da fiamme libere e fonti di calore	
Smaltire il dispositivo rispettando le normative vigenti in materia di tutela dell'ambiente e raccolta differenziata	
Confezione e Nr. Articoli	
Scatola	1 Pezzo
Garanzia	
Europea	24 mesi per difetti di fabbricazione ad esclusione delle parti consumabili
Riferimenti del Prodotto	
Fabbricante	Demarta Virginio S.r.l.
Classe di Dispositivo MDR 2017/745	Classe I
Codice ISO	06.30.18.006
Codice CND	Y063018
Codice RDM	1463825/R

Product Name

Symmetrical Prosthesis Maxima with Concaved Back.



Code	Size	French Size	European Size
06.01.4001C.0	0	65	50
06.01.4001C.1	1	70	55
06.01.4001C.2	2	75	60
06.01.4001C.3	3	80	65
06.01.4001C.4	4	85	70
06.01.4001C.5	5	90	75
06.01.4001C.6	6	95	80
06.01.4001C.7	7	100	85
06.01.4001C.8	8	105	90
06.01.4001C.9	9	110	95
06.01.4001C.10	10	115	100
06.01.4001C.11	11	120	105
06.01.4001C.12	12	125	110
06.01.4001C.13	13	130	115
06.01.4001C.14	14	135	120

Use

External breast prosthesis are used in subjects who have undergone mastectomy. Thanks to these devices the morphology of the body remains unchanged.

Description

The SYMMETRICAL PROSTHESIS MAXIMA WITH CONCAVED BACK is composed of a silicone gel wrapped in a thin polyurethane film; designed to restore an anatomical-aesthetic vision of the body. The prosthesis is hypoallergenic and aesthetically natural, the teardrop shape makes it suitable for all types of breasts.

Materials	
Prosthesis	Silicone Gel
Maintenance and Cleaning	
Cleaning is possible with a damp cloth and a mild soap	
Avoid direct contact with cosmetics (deodorant, perfumes...), please use them before putting the prosthesis on	
Storage and Disposal	
Single-user Product	
Keep away from heat, sparks, open flame or other sources of ignition	
Dispose of the device in compliance with current regulations on environmental protection and separate waste collection	
Package and Nr.Items	
Box	1 Piece
Warranty	
European	24 months for manufacturing defects with the exception of consumable parts
Product References	
Manufacturer	Demarta Virginio S.r.l.
Directive MDR 2017/745	Class I
ISO Code	06.30.18.006
CND Code	Y063018
RDM Code	1463825/R